

Illumina® Bio-Rad® SureCell™ WTA 3' Library Prep Kit for the ddSEQ™ System

Gain insight into gene expression with a sensitive, scalable, and cost-effective high-throughput next-generation sequencing platform

Highlights

- Sensitive and Unbiased Characterization of Transcriptional Signatures**
 Highly sensitive and reproducible gene detection in single cells
- Comprehensive Single-Cell RNA Sequencing Workflow**
 Fully supported workflow developed in collaboration by technology innovators
- Powerful Next-Generation Sequencing Integrated with Simple Data Analysis**
 Proven Illumina sequencing combined with streamlined, user-friendly analysis software

To deliver on the promise of single-cell biology, the Illumina Bio-Rad Single-Cell Sequencing Solution combines the highly innovative Bio-Rad Droplet Digital™ technology(3) with Illumina NGS library preparation, sequencing, and analysis technologies. This new platform provides a comprehensive, user-friendly workflow for single-cell RNA-Seq (Figure 1) that enables controlled experiments with multiple samples, treatment conditions, and time points. Built and supported in collaboration between technology leaders, the Illumina Bio-Rad Single-Cell Sequencing Solution enables transcriptome analysis of hundreds to tens of thousands of single cells in a single experiment. This scalable, robust single-cell NGS sample prep methodology enables more researchers to apply the sensitivity and precision of RNA-Seq to questions that can only be answered by interrogating individual cells.

Introduction

Complex biological systems are determined by the coordinated functions of individual cells. Conventional technologies providing bulk transcriptome data are unable to reveal the transcriptional heterogeneity that drives this complexity (1). Single-cell RNA sequencing (RNA-Seq) enables in-depth gene expression analysis, providing insight into cell function, disease progression, and therapeutic efficacy (2). However, generating thousands of single-cell next-generation sequencing (NGS) libraries in an affordable, high-throughput, and user-friendly manner remains challenging (2).

Scalable, Flexible Single-Cell Isolation

The Bio-Rad ddSEQ Single-Cell Isolator (SCI) encapsulates and partitions single cells into subnanoliter droplets on a disposable cartridge. Each cartridge can accommodate multiple samples, and multiple cartridges can be processed in parallel, allowing for isolation of hundreds to tens of thousands of cells per day. Cell lysis and cell barcoding occur inside individual droplets, enabling tracking of individual cells throughout the workflow. This droplet-based method is agnostic to mammalian cell size, providing unbiased profiling of diverse cell populations (Figure 2). The Bio-Rad ddSEQ SCI can process samples in < 5 minutes. By eliminating lengthy experimental workflows that can alter transcriptional signatures, acute transcriptional responses can be detected and tracked in time course experiments.




Fig. 1. Single-Cell RNA-Seq Workflow. The workflow integrates proven cell isolation using the Bio-Rad ddSEQ SCI, followed by library preparation using the SureCell WTA 3' Library Prep Kit with Nextera technology, Illumina sequencing, and data analysis with BaseSpace Sequence Hub and SeqGeq Analysis Software from FlowJo, LLC.




Fig. 2. Unbiased, Cell-Size Agnostic RNA-Seq. Mouse cell lines (A20, NIH 3T3) and human cell lines (HEK, BJ) were processed using the SureCell WTA 3' Library Prep Kit. Consistently high numbers of detected genes demonstrates that recovery of transcripts is not limited by cell size.

Advanced NGS Chemistry for Library Prep

A single cell suspension is loaded into a ddSEQ M Cartridge and cells are encapsulated and barcoded by the ddSEQ Single-Cell Isolator. Lysis and barcoding takes place in each droplet. Droplets are disrupted and cDNA is pooled for second strand synthesis. Libraries are generated with direct cDNA tagmentation using Nextera® technology, without the need for shearing or preamplification (Figure 3). Tagmentation is followed by 3' enrichment and sample indexing to prepare up to 24 indexed, sequencing-ready libraries with minimal hands-on time.




Fig. 3. Overview of The SureCell WTA 3' Library Prep Kit for the ddSEQ System. After cells are isolated, lysed, and barcoded, cDNA is pooled for second strand synthesis. Libraries are generated with direct cDNA tagmentation followed by 3' enrichment and sample indexing.

Reliable Sequencing with Illumina NGS Systems

Prepared single-cell libraries are loaded directly onto an Illumina MiSeq®, NextSeq®, HiSeq®, or NovaSeq™ Series System for sequencing. These sequencing systems harness industry-leading

Illumina sequencing by synthesis (SBS) chemistry, used in >90% of all NGS around the world.* This common technology provides consistency and reliability in sequencing data across all Illumina NGS platforms. The MiSeq, NextSeq, HiSeq, and NovaSeq Series of sequencing systems offer flexibility for a broad range of applications, and scalable throughput to support various study sizes.

Easy-to-use instrument control software guides researchers through sample loading and run setup with intuitive user-interfaces, and enables real-time monitoring of fully automated sequencing runs either on-instrument or remotely.

Simplified Data Analysis in BaseSpace® Sequence Hub

Single-cell sequencing data can be instantly transferred, stored, and analyzed securely in BaseSpace Sequence Hub, the Illumina cloud-based genomics computing environment. BaseSpace Sequence Hub provides a large collection of BaseSpace Apps. Commercial and open-source tools support a range of common data analysis needs such as alignment, variant calling, and more. These Apps feature intuitive push-button user interfaces designed to be used without the need for bioinformatics expertise.

The SureCell RNA Single-Cell App was designed to support data analysis for the Illumina Bio-Rad Single-Cell Sequencing Solution. The SureCell RNA Single-Cell App enables streamlined data analysis for up to 96 samples across multiple sequencing runs and includes:

- Convenient sequencing quality control (QC) metrics
- Easy assignment of unique transcripts to single cells
- Exportable and downloadable gene expression matrices and reports
- Various analysis options for identification of subpopulations and differentially expressed genes

Advanced Data Visualization with SeqGeq Software



SeqGeq Software is a desktop application for advanced data analysis, exploration, and visualization of single-cell gene expression data developed by FlowJo, LLC for the Illumina Bio-Rad Single-Cell Sequencing Solution. SeqGeq Software features powerful data reduction and population identification tools. Direct integration with BaseSpace Sequence Hub enables visualization and analysis of expression data with statistic color-mapping of individual cells, summary heatmaps, and drag-and-drop report editors.

High-Quality Data

To demonstrate the high-quality single-cell RNA-Seq data achieved with the Illumina Bio-Rad Single-Cell Sequencing Solution, *in vitro* experiments were performed with cells mixed from different species. Human embryonic kidney 293 (HEK 293) cells and NIH 3T3 mouse

*Data calculations on file. Illumina, Inc. 2015.

embryonic fibroblasts were mixed at a 1:1 ratio (unless otherwise noted) and loaded across four sample chambers of a single Bio-Rad ddSEQ M Cartridge. The ddSEQ SCI encapsulated and barcoded 1,384 single cells. Barcoded transcripts were processed for single-cell sequencing using the Illumina Bio-Rad SureCell WTA 3' Library Prep Kit for the ddSEQ System, and sequenced on the Illumina NextSeq 550 System. Sequencing results were analyzed using the Illumina SureCell RNA Single-Cell App.

Confident Gene Detection in a Heterogeneous Population of Cells

Data analysis and plotting of human (hg19) and mouse (mm10) unique transcripts detected in the mixed species samples shows a low percentage of cell barcodes (5.8%) containing transcripts mapping to both species, which represent cell doublets (Figure 4A). This demonstrates efficient single-cell isolation at high purity (99%). By graphing the cumulative fraction of gene transcripts assigned to cell barcodes, the inflection point (red line) is used to determine the number of barcoded cells detected (Figure 4B). This indicates that a high fraction of transcripts are assigned to single cells.




Fig. 4. Confident Single-Cell Transcript Identification. **A.** BaseSpace generated plots of the number of unique molecular identifiers (UMI), ie transcripts, assigned to the mouse/mm10 (red) and human/hg19 (blue) genomes for each cell barcode. Unique transcripts mapping to both mouse and human (purple) represent cell doublets. **B.** Cumulated fraction of unique transcripts assigned to cell barcodes (linear scale). The inflection point (red line) determines the number of barcoded cells detected in the sequencing run.

Clear Identification of Cell Subpopulations in Mixed Samples

An *in vitro* species mixing experiment was performed in which mouse NIH 3T3 cells were spiked with a small number of human HEK 293 cells. With t-distributed stochastic neighbor embedding (t-SNE) analysis, human cells from this mixture were identified as a distinct cluster representing 7% of the total cell population (Figure 5A). Cells color coded by gene expression of human *RPL 13* confirm the identity of the subpopulation as human cells (Figure 5B).




Fig. 5. Confident Cell Subpopulation Identification. **A.** Analysis in the SureCell RNA Single-Cell App using the t-SNE algorithm of a mixture of NIH 3T3 and HEK 293 cells identifies a distinct subpopulation of cells. **B.** Cells color coded by gene expression of hg19 *RPL 13* confirms the identity of the subpopulation as human.

Analysis of Single Cells Reveals Heterogeneity of Cell Cycle Phases

An *in vitro* species mixing experiment was performed in which HEK 293 cells were mixed at a 1:1 ratio with NIH 3T3 cells. Principal component analysis (PCA) and color-mapped gene expression plot of human *RPL 13* identifies human cells (Figure 6A). t-SNE analysis of the subset of HEK 293 cells reveals distinct clusters of cells, representing differences in gene expression profiles (Figure 6B). Deeper analysis of markers associated with phases of the cell cycle enables accurate cell cycle analysis of single cells. This analysis is based on unique transcript counts of relevant genes, normalized by total count for each cell and visualized as a summary heatmap (Figure 6C).




Fig. 6. Deep Single-Cell Analysis in SeqGeq Software Reveals Homogeneity in HEK293 Cells. PCA of 1:1 ratio mixture of HEK293 and NIH3T3 cells, color-mapped by expression of human *RPL13* gene (A, mouse cells in gray). t-SNE analysis of the human cell subset where distinct clusters represent differences in gene expression profiles (B). Single cell heatmap reveals cell cycle state based on unique transcript counts of genes (C). Expression is centered by the median and scaled by the median absolute deviation for each cell cycle. Analysis was performed using BaseSpace Sequence Hub and SeqGeq Software.




Fig. 7. Sensitivity and Reproducibility of Gene Detection Across Varied Cell Lines. **A.** The median number of genes detected per cell are plotted at sequencing depths from 25,000–200,000 reads per cell for two replicates of human HEK 293 cells. **B.** The median number of genes detected per cell for two replicates of mouse NIH 3T3 cells. **C.** A linear regression of gene expression is plotted for genes with ≥ 50 counts, summed across all HEK 293 cells in two replicate samples.

Highly Sensitive and Reproducible Results

Replicate samples of HEK 293 cells were processed and sequenced on a NextSeq 550 System as described. Sequencing reads were subsampled to various reads per cells ranging from 25,000–200,000 reads. The median number of reads detected per cell at each sequencing depth shows highly sensitive and reproducible gene detection across replicates of different cell lines (Figure 7A and 7B). Total gene counts for each human gene were summed across all HEK

293 cells. A linear regression fit of the summed gene counts between two replicates processed on a single ddSEQ M Cartridge further demonstrates highly reproducible results (Figure 7C).

Summary

To advance understanding of the transcriptional heterogeneity that drives complex biological systems, researchers need a scalable, high-throughput, and user-friendly method for generating thousands of single-cell NGS libraries. The Illumina Bio-Rad Single-Cell Sequencing Solution is a comprehensive workflow developed in collaboration by industry experts in droplet-based cell isolation and NGS technologies. It reveals new types of single-cell information by facilitating analysis of multiple samples in parallel, under multiple treatment conditions, and at multiple time-points. The simple yet powerful data analysis options in the SureCell RNA Single-Cell App, combined with data reduction and population identification tools in SeqGeq Software, can resolve heterogeneous cell populations and identify subpopulations of interest using gene expression profiles and data visualization tools. Analysis of cell cycle markers allows for cell cycle analysis of individual cells in complex tissues. The Illumina Bio-Rad Single-Cell Sequencing Solution enables highly sensitive and reproducible interrogation of single-cell transcriptomes from hundreds to tens of thousands of single cells in a single experiment.

Ordering Information

Product	Catalog No.
SureCell WTA 3' Library Prep Kit (2 cartridge kit)	20014279
SureCell WTA 3' Library Prep Kit (6 cartridge kit)	20014280

Learn More

To learn more about the Illumina Bio-Rad Single-Cell Sequencing Solution, visit:

www.illumina.com/surecell

www.bio-rad.com/ddSEQ

www.flowjo.com/seqgeq

References

1. Kolodziejczyk AA, Kim JK, Svensson V, Marioni JC, Teichmann SA. The technology and biology of single-cell RNA sequencing. *Mol Cell*. 2015;58(4):610–620.
2. Shapiro E, Biezuner T, Linnarsson S. Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet*. 2013;14(9):618–630.
3. Hindson BJ, Ness, KD, Masquelier DA, et al. High-throughput droplet digital PCR system for absolute quantitation of DNA copy number. *Anal Chem*. 2011;83(22):8604–8610.

For Research Use Only. Not for use in diagnostic procedures.