




Small Genome Workflow on the MiniSeq™ System



This checklist supports the following workflow choices.

Sequencing Instrument:	MiniSeq System
Setup Option	Local Run Manager
Library Preparation Kit:	Nextera XT DNA Sample Prep
Indexing:	Dual Indexing
Reagent Kits:	MiniSeq Kit
Analysis Workflow:	Resequencing
Analysis Software:	Local Run Manager

Small Genome Checklist

Set Run Parameters

- 1 Log in to Local Run Manager.
- 2 Click **Create Run**, and select **Resequencing**.
- 3 Enter a run name that identifies the run.
- 4 [Optional] Enter a run description.
- 5 From the Library Kit drop-down list, select Nextera XT or Nextera XT V2.
- 6 Click **2** to specify a dual-indexed run.
- 7 Specify a read type: **Single Read** or **Paired End**.
- 8 Enter the number of cycles for the run.
- 9 Select an alignment method.
 - ▶ **BWA-MEM**—(Default) Optimized for Illumina sequencing data and reads ≥ 70 bp.
 - ▶ **BWA-Backtrack Legacy**—Use with legacy data or reads < 70 bp.
- 10 Select a variant calling method.
 - ▶ **Starling**—(Default) Calls SNPs and small indels, and summarizes depth and probabilities for every site in the genome.
 - ▶ **GATK**—Calls raw variants for each sample, analyzes variants against known variants, and then calculates a false discovery rate for each variant.
- 11 Enable or disable the Export gVCF, Flag PCR Duplicates, and Indel Realignment settings.
- 12 Enter a unique sample ID.
- 13 [Optional] Enter a sample description.
- 14 Select an Index 1 adapter.
- 15 Select an Index 2 adapter.
- 16 Select a reference genome.
- 17 Click **Save Run**.

Tagment Genomic DNA

- 1 Add the following to a new PCR plate. Pipette to mix.

Item	Volume (μ l)
TD	10
Normalized gDNA	5
- 2 Add 5 μ l ATM. Pipette to mix.
- 3 Centrifuge at $280 \times g$ at 20°C for 1 minute.
- 4 Place on the thermal cycler and run the tagmentation program.
- 5 Add 5 μ l NT. Pipette to mix.
- 6 Centrifuge at $280 \times g$ at 20°C for 1 minute.
- 7 Incubate at room temperature for 5 minutes.
- 8 [Optional] Run 1 μ l sample on a High Sensitivity DNA chip.

Amplify Libraries

- 1 [24 libraries] Arrange the index primers as follows.
 - ▶ Arrange Index 1 (i7) adapters in columns 1–6.
 - ▶ Arrange Index 2 (i5) adapter in rows A–D.
- 2 [96 libraries] Arrange the index primers as follows.
 - ▶ Arrange Index 1 (i7) adapters in columns 1–12.
 - ▶ Arrange Index 2 (i5) adapter in rows A–H.
- 3 Add 5 µl of each Index 1 (i7) adapter down each column. Replace the cap on each i7 adapter tube with a new orange cap.
- 4 Add 5 µl of each Index 2 (i5) adapter across each row. Replace the cap on each i5 adapter tube with a new white cap.
- 5 Add 15 µl NPM. Pipette to mix.
- 6 Centrifuge at 280 × g at 20°C for 1 minute.
- 7 Place on the thermal cycler and run the PCR program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

Clean Up Libraries

- 1 Centrifuge at 280 × g at 20°C for 1 minute.
- 2 Transfer 50 µl PCR product.
- 3 Add 30 µl AMPure XP beads.
- 4 Shake at 1800 rpm for 2 minutes.
- 5 Incubate at room temperature for 5 minutes.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Wash 2 times with 200 µl 80% EtOH.
- 9 Using a 20 µl pipette, remove residual 80% EtOH.
- 10 Air-dry on the magnetic stand for 15 minutes.
- 11 Remove from the magnetic stand.
- 12 Add 52.5 µl RSB.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 50 µl supernatant.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Check Libraries

- 1 [Optional] Run 1 µl of undiluted library on an Agilent Technology 2100 Bioanalyzer using a High Sensitivity DNA chip.

Normalize Libraries

- 1 Transfer 20 µl supernatant.
- 2 [96 samples] Add 4.4 ml LNA1 to a new 15 ml conical tube.
- 3 Thoroughly resuspend LNB1. Pipette to mix.
- 4 Transfer 800 µl LNB1 to the tube. Invert to mix.
- 5 Add 45 µl combined LNA1/LNB1.
- 6 Shake at 1800 rpm for 30 minutes.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 45 µl LNW1.
- 10 Add 30 µl 0.1 N NaOH.
- 11 Shake at 1800 rpm for 5 minutes.
- 12 During the 5 minute elution, label a new 96-well PCR plate SGP.
- 13 Add 30 µl LNS1 to the SGP plate. Set aside.
- 14 After the 5 minute elution, make sure that all samples are resuspended. Pipette to mix.
- 15 Shake at 1800 rpm for 5 minutes.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer the supernatant from the midi plate to the SGP plate.
- 18 Centrifuge at 1000 × g for 1 minute.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Pool Libraries

- 1 Centrifuge at 1000 × g at 20°C for 1 minute.
- 2 Transfer 5 µl from the SGP plate to a new PCR 8-tube strip.
- 3 Label a new Eppendorf tube PAL.
- 4 Transfer the contents of the PCR 8-tube strip to the PAL tube. Invert to mix.
- 5 Dilute pooled libraries to the loading concentration for the sequencing instrument you are using. See the denature and dilute libraries guide for your instrument.
- 6 Store unused pooled libraries in the PAL tube and SGP plate at -25°C to -15°C for up to 7 days.

Prepare Consumables

- 1 Remove the reagent cartridge from -25°C to -15°C storage.
- 2 Thaw reagents in a room temperature water bath for 90 minutes.
- 3 Invert the cartridge 5 times to mix reagents.
- 4 Gently tap on the bench to reduce air bubbles.
- 5 Remove a new flow cell package from 2°C to 8°C storage.
- 6 Set the unopened flow cell package aside at room temperature for 30 minutes.
- 7 Remove the flow cell from the foil package and flow cell container.
- 8 Clean the glass surface of the flow cell with a lint-free alcohol wipe.
- 9 Dry with a lint-free lens cleaning tissue.

Denature, Dilute, and Load Libraries

- 1 Thaw the Hybridization Buffer at room temperature.
- 2 Vortex briefly before use.
- 3 Preheat the incubator to 98°C.
- 4 Combine the 5 µl pooled libraries and 995 µl prechilled Hybridization Buffer in a microcentrifuge tube.
- 5 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 6 Transfer 250 µl diluted library to a new microcentrifuge tube.
- 7 Add 250 µl prechilled Hybridization Buffer.
- 8 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 9 Place the tube on the preheated incubator for 2 minutes.
- 10 Immediately cool on ice.
- 11 Leave on ice for 5 minutes.
- 12 [Optional] Denature and dilute a PhiX control to 1.8 pM and a 1% spike-in to the final library.
- 13 Clean the foil seal covering reservoir #16 using a low-lint tissue.
- 14 Pierce the seal with a clean 1 ml pipette tip.
- 15 Add 500 µl prepared 1.8 pM libraries into reservoir #16.

Perform a Sequencing Run

- 1 From the Home screen, select **Sequence**.
- 2 Enter your user name and password.
- 3 Select **Next**.
- 4 Select a run name from the list of available runs.
- 5 Select **Next**.
- 6 Open the flow cell compartment door.
- 7 Press the release button to the right of the flow cell latch.
- 8 Place the flow cell on the flow cell stage over the alignment pins.
- 9 Close the flow cell latch to secure the flow cell.
- 10 Close the flow cell compartment door.
- 11 Open the reagent compartment door.
- 12 Slide the reagent cartridge into the reagent compartment until the cartridge stops.
- 13 Remove the spent reagents bottle from the compartment.
- 14 Discard the contents and slide the empty spent reagents bottle into the compartment.
- 15 Close the compartment door and select **Next**.
- 16 Confirm run parameters.
- 17 Select **Next**.
- 18 When the automated check is complete, select **Start**.
- 19 Monitor run progress, intensities, and quality scores as metrics appear on the screen.

View Analysis Results

- 1 From the Local Run Manager dashboard, click the run name.
- 2 From the Run Overview tab, review the sequencing run metrics.
- 3 [Optional] Click the **Copy to Clipboard**  icon for access to the output run folder.
- 4 Click the Sequencing Information tab to review run parameters and consumables information.
- 5 Click the Samples and Results tab to view the analysis report.
- 6 [Optional] Click the **Copy to Clipboard**  icon for access to the Analysis folder.