

Illumina Free Adapter Blocking Reagent

Reference Guide

Introduction	3
Sample Treatment Workflow	3
Sample Treatment	4
Technical Assistance	6



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2018 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

Protocol

Introduction	3
Sample Treatment Workflow	3
Sample Treatment	4

Introduction

This protocol provides information for performing a free adapter blocking treatment on libraries to block excess free adapter, minimizing potential index hopping levels and enhancing data quality.

Use this method on unpooled or pooled libraries with any NGS library prep method before the denaturation step and subsequent cluster generation using Illumina sequencers.

Sample Treatment Workflow

The following diagram illustrates the treatment workflow.

Figure 1 Sample Treatment Workflow



Sample Treatment

The following consumables and equipment are required to perform the protocol. Part numbers are provided where applicable.

Consumables and Equipment

Item	Part Number	Supplied By
Illumina Free Adapter Blocking Reagent (FAB) 12 reactions	20024144	Illumina
Illumina Free Adapter Blocking Reagent (FAB) 48 reactions	20024145	Illumina
Agencourt Ampure XP beads	A63880	Beckman Coulter Life Sciences
Library sample pool – 30 µl at concentration to be used for denaturation during clustering		User
Freshly prepared 80% ethanol (EtOH)		User
Ethanol 200 proof, for molecular biology (absolute)	E7023	Sigma-Aldrich
Magnetic Stand-96	AM10027	Thermo Fisher Scientific
Vortexer		General lab supplier
96-well thermal cycler (with heated lid)		General lab supplier

About Reagents

- ▶ Vortex Ampure XP before each use.
- ▶ Vortex Ampure XP frequently to make sure that beads are evenly distributed.
- ▶ Aspirate and dispense Ampure XP slowly due to the viscosity of the solution.

Preparation

- 1 Prepare the following consumables.

Item	Storage	Instructions
FAB	-25°C to -15°C	Thaw at room temperature, and then place on ice. Return to storage after use.
RSB	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.
Ampure XP beads	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.

- 2 Save the FAB program on the thermal cycler:
 - ▶ Choose the preheat lid option and set to 100°C.
 - ▶ 38°C for 20 minutes
 - ▶ 60°C for 20 minutes
 - ▶ Hold at 4 °C
 Each tube contains 60 µl.

Procedure

Sample Treatment

- 1 Mix FAB reagent by inversion.
- 2 Centrifuge FAB at $600 \times g$ for 5 seconds.
- 3 Add 30 μ l of library sample pool to PCR tube.
- 4 Add 30 μ l of FAB to each PCR tube and then mix thoroughly by pipetting up and down.
- 5 Centrifuge briefly to make sure all contents are on the bottom of the tube.
- 6 Incubate by placing on the thermal cycler and running the FAB program.

Clean Up Pooled Sample

- 1 Vortex AMPure XP until well-dispersed.
- 2 Add 60 μ l AMPure XP to each sample treatment tube and mix thoroughly by pipetting up and down.
- 3 Incubate at room temperature for 5 minutes.
- 4 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 5 Remove and discard all supernatant from each tube.
- 6 Wash 2 times as follows:
 - ▶ Add 200 μ l freshly prepared 80% EtOH to each tube.
 - ▶ Incubate on the magnetic stand for 30 seconds.
 - ▶ Remove and discard all supernatant from each tube.
- 7 Use a 20 μ l pipette to remove residual EtOH from each tube.
- 8 Air-dry on the magnetic stand for 5 minutes.
- 9 Add 22.5 μ l RSB to each tube.
- 10 Remove from the magnetic stand and then mix thoroughly by pipetting up and down.
- 11 Incubate at room temperature for 2 minutes.
- 12 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 13 Transfer 20 μ l supernatant to a new tube.
- 14 Quantify libraries if necessary and proceed onto standard clustering for your Illumina sequencing instrument starting with the NaOH denaturation step.
- 15 Store at -25°C to -15°C if not clustering immediately.

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com
Email: techsupport@illumina.com

Illumina Customer Support Telephone Numbers

Region	Toll Free	Regional
North America	+1.800.809.4566	
Australia	+1.800.775.688	
Austria	+43 800006249	+43 19286540
Belgium	+32 80077160	+32 34002973
China	400.066.5835	
Denmark	+45 80820183	+45 89871156
Finland	+358 800918363	+358 974790110
France	+33 805102193	+33 170770446
Germany	+49 8001014940	+49 8938035677
Hong Kong	800960230	
Ireland	+353 1800936608	+353 016950506
Italy	+39 800985513	+39 236003759
Japan	0800.111.5011	
Netherlands	+31 8000222493	+31 207132960
New Zealand	0800.451.650	
Norway	+47 800 16836	+47 21939693
Singapore	+1.800.579.2745	
Spain	+34 911899417	+34 800300143
Sweden	+46 850619671	+46 200883979
Switzerland	+41 565800000	+41 800200442
Taiwan	00806651752	
United Kingdom	+44 8000126019	+44 2073057197
Other countries	+44.1799.534000	

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download in PDF from the Illumina website. Go to support.illumina.com, select a product, then select **Documentation & Literature**.



Illumina
5200 Illumina Way
San Diego, California 92122 U.S.A.
+1.800.809.ILMN (4566)
+1.858.202.4566 (outside North America)
techsupport@illumina.com
www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2018 Illumina, Inc. All rights reserved.

illumina[®]