

Fragment DNA

- 1 Normalize gDNA to 10 ng/μl. Pipette or vortex to mix.
- 2 Centrifuge briefly.
- 3 Transfer 50 μl DNA to Covaris tubes.
- 4 Centrifuge at 280 × g for 5 seconds.
- 5 Fragment the DNA.
- 6 Centrifuge at 280 × g for 5 seconds.
- 7 Transfer 50 μl sample volume to a new eight-tube strip.
- 8 [Recommended] Run 1 μl using a High Sensitivity DNA chip.

SAFE STOPPING POINT

If you are stopping, cap the tubes and store at 2°C to 8°C overnight.

Clean Up Fragmented DNA

- 1 Add 80 μl SPB. Pipette or vortex to mix.
- 2 Incubate at room temperature for 3 minutes.
- 3 Place on a magnetic stand until liquid is clear.
- 4 Remove and discard all supernatant.
- 5 Wash two times with 200 μl 80% EtOH.
- 6 Centrifuge briefly, and then place on a magnetic stand.
- 7 Remove and discard residual EtOH.
- 8 Air-dry for 5 minutes.
- 9 Add 60 μl RSB.
- 10 Remove from the magnetic stand, and then pipette or vortex to mix.
- 11 Incubate at room temperature for 1 minute.

Repair Ends

- 1 Invert ERP3 to mix, and then tap.
- 2 Add 40 μl ERP3. Pipette or vortex to mix.
- 3 Place on the thermal cycler and run the ERP program. Remove when samples reach 4°C.
- 4 Briefly centrifuge at 280 × g.
- 5 Add 120 μl SPM.
- 6 Pipette to mix.
- 7 Incubate at room temperature for 3 minutes.
- 8 Place on a magnetic stand until liquid is clear. Keep on the magnetic stand until step 11.
- 9 Remove and discard supernatant.
- 10 Wash two times with 200 μl 80% EtOH.
- 11 Centrifuge briefly.
- 12 Place on a magnetic stand. Remove and discard residual EtOH.
- 13 Air-dry for 5 minutes.
- 14 Remove from the magnetic stand.
- 15 Add 17.5 μl RSB. Pipette or vortex to mix.
- 16 Incubate at room temperature for 1 minute.

Adenylylate 3' Ends

- 1 Add 12.5 µl ATL2.
- 2 Pipette or vortex to mix.
- 3 Place on the thermal cycler and run the ATAIL70 program.

Ligate Adapters

- 1 Briefly centrifuge the DNA adapter tubes.
- 2 Remove LIG2 from storage.
- 3 Add the following to the 8-tube strip containing sample.
 - ▶ RSB (1.5 µl)
 - ▶ LIG2 (2.5 µl)
 - ▶ DNA adapters (4 µl)
- 4 Pipette or vortex to mix, and then centrifuge briefly.
- 5 Place on the thermal cycler and run the LIG program. Remove when samples reach 4°C.
- 6 Centrifuge briefly.
- 7 Vortex STL to mix.
- 8 Add 5 µl STL. Pipette or vortex to mix.
- 9 Add 43 µl SPM. Pipette or vortex to mix.
- 10 Incubate at room temperature for 3 minutes.
- 11 Place on a magnetic stand until liquid is clear. Keep on the magnetic stand until step .
- 12 Remove and discard supernatant.
- 13 Wash two times with 200 µl 80% EtOH.
- 14 Centrifuge briefly, and then place on a magnetic stand.
- 15 Remove and discard residual EtOH.
- 16 Air-dry for 5 minutes.
- 17 Elute and combine sets of four samples with 40 µl RSB.

Hybridize Probes

- 1 Add the following reagents in the order listed to the tube.
 - ▶ Blocker (10 µl)
 - ▶ EPIC Oligos (10 µl)
 - ▶ SPM (150 µl)
- 2 Pipette carefully to mix, and make sure that the mixture is homogenous.
- 3 Incubate at room temperature for 10 minutes.
- 4 Place on a magnetic stand until liquid is clear.
- 5 Remove and discard supernatant.
- 6 Wash two times with 200 µl 80% EtOH.
- 7 Centrifuge briefly.
- 8 Place on a magnetic stand. Remove and discard residual EtOH.
- 9 Air-dry for 5 minutes.
- 10 Add 7.7 µl CT4. Pipette or vortex to mix.
- 11 Incubate at room temperature for 2 minutes.
- 12 Place on a magnetic stand until liquid is clear.
- 13 Add 2.5 µl EHB2 to a new tube.
- 14 Transfer 7.5 µl to the eight-tube strip containing 2.5 µl EHB2.
- 15 Pipette or vortex to mix.
- 16 Place on the thermal cycler and run the MC HYB1 program.
- 17 Keep at the 58°C holding temperature for 35 minutes–2 hours.

Capture Hybridized Probes

- 1 Centrifuge tubes at 280 × g for 1 minute.
- 2 Vortex SMB to mix.
- 3 Add 250 µl SMB to a new tube.
- 4 Transfer each pool to the tubes of SMB.
- 5 Rinse out hybridization tube with SMB.
- 6 Vortex to mix.
- 7 Incubate at room temperature for 25 minutes.
- 8 Centrifuge briefly.
- 9 Place on a magnetic stand until liquid is clear.
- 10 Remove and discard all supernatant.
- 11 Remove from the magnetic stand.
- 12 Vortex EWS.
- 13 Add 200 µl EWS.
- 14 Pulse vortex to mix until SMB pellet is resuspended.
- 15 Place on the 50°C heat block for 20 minutes.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Remove and discard all supernatant.
- 18 Remove from the magnetic stand.
- 19 Repeat steps 13–18.
- 20 Centrifuge briefly.
- 21 Place on a magnetic stand and remove the remaining supernatant.
- 22 Create elution premix in a 1.7 ml microcentrifuge tube.
 - ▶ ET1 (31.4 µl)
 - ▶ HP3 (1.7 µl)
- 23 Add 30 µl elution premix. Pipette or vortex to mix.
- 24 Incubate at room temperature for 5 minutes, and then centrifuge briefly.
- 25 Place on a magnetic stand until liquid is clear.
- 26 Add 5 µl ET2 to a new eight-tube strip.

- 27 Transfer 29 µl supernatant to ET2. Pipette or vortex to mix.

Perform Second Hybridization

- 1 Add 6 µl water, 50 µl CT3, and 10 µl EPIC Oligos to the tube from the previous step.
- 2 Pipette or vortex to mix, and then centrifuge briefly.
- 3 Place on the thermal cycler and run the MC HYB2 program.

SAFE STOPPING POINT

Keep at the 58°C holding temperature for at least 14.5 hours.

Perform Second Capture

- 1 Centrifuge tubes at 280 × g for 1 minute.
- 2 Vortex SMB to mix.
- 3 Add 250 µl SMB to a new tube.
- 4 Transfer each pool to the tubes.
- 5 Vortex to mix.
- 6 Incubate at room temperature for 25 minutes.
- 7 Centrifuge briefly.
- 8 Place on a magnetic stand until liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Remove from the magnetic stand.
- 11 Vortex EWS.
- 12 Add 200 µl EWS.
- 13 Pulse vortex to mix.
- 14 Place on the 50°C heat block for 30 minutes.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Remove and discard all supernatant.
- 17 Remove from the magnetic stand.
- 18 Repeat steps 12–17.
- 19 Centrifuge briefly.
- 20 Place on a magnetic stand and remove the remaining supernatant.
- 21 Create elution premix in a 1.7 ml microcentrifuge tube.
 - ▶ ET1 (18.8 µl)
 - ▶ HP3 (1 µl)
- 22 Add 18 µl elution premix.
- 23 Pipette or vortex to mix.
- 24 Incubate at room temperature for 5 minutes, and then centrifuge briefly.
- 25 Place on a magnetic stand until liquid is clear.
- 26 Add 2.9 µl ET2 to a new eight-tube strip.
- 27 Transfer 17.1 µl supernatant to the tube.
- 28 Pipette or vortex to mix. Proceed *immediately*.

Bisulfite Conversion

- 1 Add 130 µl Lightning Conversion Reagent to sample.
- 2 Pipette or vortex to mix, and then centrifuge briefly.
- 3 Place on the thermal cycler and run the BSF CON program.
- 4 Add 600 µl M-Binding Buffer and 10 µl MagBinding Beads to a new tube.
- 5 Transfer entire sample to the tube of M-Binding Buffer and MagBinding Beads, and then vortex to mix.
- 6 Incubate at room temperature for 5 minutes.
- 7 Centrifuge briefly.
- 8 Place on a magnetic stand until liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Wash one time with 400 µl M-Wash Buffer.
- 11 Add 200 µl L-Desulphonation Buffer, and then vortex to mix.
- 12 Incubate at room temperature for 15 minutes.
- 13 Place on a magnetic stand until liquid is clear.
- 14 Remove and discard all supernatant.
- 15 Wash two times with 400 µl M-Wash Buffer.
- 16 Centrifuge briefly, and then place on a magnetic stand.
- 17 Remove and discard supernatant.
- 18 Place on the 50°C heat block until dry.
- 19 Add 23 µl RSB. Pipette or vortex thoroughly to mix.
- 20 Place on the 50°C heat block for 4 minutes.
- 21 Place on a magnetic stand until liquid is clear.
- 22 Transfer 20 µl DNA to a new eight-tube strip.

Amplify Enriched Library

- 1 Add the following reagents to 20 µl DNA from the previous step.
 - ▶ PPC (5 µl)
 - ▶ Kapa HiFi Uracil+ (25 µl)
- 2 Pipette or vortex to mix.
- 3 Place on the thermal cycler and run the AMP MC program.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

Clean Up Amplified Enriched
Library

- 1 Add 50 µl SPB to the tube, and then pipette or vortex to mix.
- 2 Incubate at room temperature for 3 minutes.
- 3 Place on a magnetic stand until liquid is clear.
- 4 Remove and discard supernatant.
- 5 Wash two times with 200 µl 80% EtOH.
- 6 Centrifuge briefly.
- 7 Place on a magnetic stand and remove residual EtOH.
- 8 Air-dry until dry.
- 9 Add 20 µl RSB, and then pipette or vortex to mix.
- 10 Incubate at room temperature for 1 minute.
- 11 Place on a magnetic stand until liquid is clear.
- 12 Transfer 19 µl supernatant to a new tube.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 6 months.

Check Enriched Libraries

- 1 Quantify the postenriched library using the Qubit dsDNA HS Assay Kit.
 - ▶ Use 2 µl as the loading volume.
 - ▶ Use the dsDNA and high sensitivity settings.
 - ▶ Record STD1 and STD2 readings.
 - ▶ Measure the library concentration.
- 2 Run 1 µl using a High Sensitivity DNA chip.

Acronyms

Acronym	Definition
ATL2	A-Tailing Mix
BLR	Blocker
CT3	Capture Target Buffer 3
CT4	Capture Target Buffer 4
EDTA	Ethylenediaminetetraacetic Acid
EPIC	Epigenetic Oligo Pool
ERP2	End Repair Mix 2
ET1	Elute Target Buffer 1
ET2	Elute Target Buffer 2
EWS	Enrichment Wash Solution
HP3	2 N NaOH
LIG2	Ligation Mix 2
PPC	PCR Primer Cocktail
RSB	Resuspension Buffer
SMB	Streptavidin Magnetic Beads
SPB	Sample Purification Beads
SPM	Sample Purification Mix
STL	Stop Ligation Buffer