cBot 2 System Safety and Compliance Guide

For Research Use Only. Not for use in diagnostic procedures.

This guide provides important safety information pertaining to the installation, servicing, and operation of the Illumina® cBot™ 2 system. Read this document before performing any procedures on the cBot 2.

The country of origin and date of manufacture of the system are printed on the instrument label.

Safety Considerations and Markings

This section identifies potential hazards associated with installing, servicing, and operating the instrument. Do not operate or interact with the instrument in a manner that exposes you to any of these dangers.

All of the hazards described herein can be avoided by following the standard operating procedures included in the *cBot 2 System Guide* (document # 15065681).

General Safety Warnings

Make sure that all personnel are trained in the correct operation of the instrument and any potential safety considerations.

Follow all operating instructions when working in areas marked with this label to minimize risk to personnel or the instrument.

Laser Safety Warning

The barcode scanner on the front of the cBot is a Class 2 laser product.

Do not stare into the visible-light beam of the barcode scanner.

The instrument is powered by 100–240 VAC operating at 50–60 Hz. Power consumption is 500 W. Hazardous voltage sources are located behind the right side panel, but can be accessible if other panels are removed. Some voltage is present on the instrument even when the instrument is turned off. Operate the instrument with all panels intact to avoid electrical shock.

Protective Earth

The instrument has a connection to protective earth through the enclosure. The safety ground on the power cord returns protective earth to a safe reference. The protective earth connection on the power cord must be in good working condition when using this device.

Fuses

The instrument contains no user-replaceable fuses. The instrument uses only manufacturer-recommended fuses. These fuses are size 5x20 and are rated for 8 Amps, 250 VAC, Slo-Blo.

Hot Surface Safety Warning

Do not operate the instrument with any of the panels removed.

Do not touch the aluminum thermal block on the thermal stage. The heater used in this area is normally controlled between ambient room temperature (22°C) and 95°C. Exposure to temperatures at the upper end of this range can result in burns. Operate the instrument only when the lid closed.

Pinch Hazard Warning

Do not open the cBot lid when the instrument is running. Mechanical moving parts in the reagent compartment can injure fingers and hands. Keep fingers and hands clear of the area around the reagent stage, and do not reach behind the sipper comb.

Environmental Considerations

Element	Specification
Temperature	Maintain a lab temperature of 19°C to 25°C (22°C ±3°C). This temperature is the operating temperature of the instrument.
Humidity	Maintain a noncondensing relative humidity between 20–80%.
Elevation	Locate the instrument at an altitude below 2000 meters (6500 feet).
Air Quality	Operate the instrument in a Pollution Degree II environment or better. A Pollution Degree II environment is defined as an environment that normally includes only nonconductive pollutants.
Ventilation	Maximum thermal output is approximately 1700 BTU/h (500 W).

Noise Output

Noise output is 65 dB when standing 1 m (3.3 ft) from the front of the instrument.

Product Certifications and Compliance

The cBot 2 is certified to the following electrical standards:

- ▶ Conforms to UL STD 61010-1
- ▶ Certified to CSA STD C22.2 No 61010-1
- ▶ Low Voltage Directive 2006/95/EC
- ▶ EMC Directive 2004/108/EC
- ▶ IEC/EN 61010-1
- ▶ IEC/EN 61326-1
- ▶ IEC/EN 60825-1

Compliance and Regulatory Markings

The instrument is labeled with the following compliance and regulatory markings.

This label assures that the product is compliant with safety standards and has been tested and certified by a third-party organization.

This label assures that the product meets the essential requirements of all relevant EU directives.

This label assures that the product complies with the Environmental Protection User Period - 10 years.

Restriction of Hazardous Substances (RoHS)

This label indicates that the instrument meets the WEEE Directive for waste.

Return the instrument to Illumina for disposal.

Product Compliance and Regulatory Statements

Laser Safety Labels

The cBot 2 is labeled with the following labels for safety.

Figure 1 Class 2 Laser Warning


Figure 2 Class 1 Laser Warning

Laser Safety Labels (French Translations)

Figure 3 Class 2 Laser Warning (French)


Figure 4 Class 1 Laser Warning (French)

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Notes

Revision History

Document	Date	Description of Change
Material # 20000254 Document # 15065643 v01	January 2016	Added translations for French, German, Italian, Korean, Russian, and Spanish.
Part # 15065643	May 2015	Initial release.

Copyright and Trademarks

© 2016 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.