

Clarity LIMS v6.0.1

Release Notes

February 9, 2023

Introduction

These Release Notes describe the key changes made to software components for Clarity LIMS since version 6.0. When upgrading from a version prior to v6.0, review the release notes for v6.0 for a list of features and bug fixes introduced in that version. This is a required patch. Clarity LIMS v6.0.1 is deployable in both hosted and on-premise environments.

I. Clarity v6.0.1

NEW FEATURES:

- None

DEFECT REPAIRS:

- Fixed an issue that impacted the link between the event log and audit trail.
- Fixed an issue that resulted in the unexpected addition of controls to workflows containing related samples.

KNOWN ISSUES:

- None