

AutoConvert Software 2.0.1

Release Notes

Introduction

These Release Notes detail the key changes to the AutoConvert Software since the release of AutoConvert Software v2.0.0. This is a required software update for genotyping customers currently using AutoConvert Software v2.0.0, as this release includes important features and bug fixes outlined below.

IMPROVEMENTS:

- Improved gender estimation algorithm to improve accuracy in situations where the heterozygosity of SNPs on the X chromosome is low due to low average minor allele frequencies (MAF) or when the sample has a substantial loss of heterozygosity (LOH).

KNOWN ISSUES:

- None