



Beeline v2.0.3

Release Notes

FOR RESEARCH USE ONLY

© 2017 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, GAIIX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, Nextera, NextBio, NextSeq, NovaSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.



Introduction

These Release Notes detail the key changes to software components for Beeline v2.0.3 since the v2.0.2 release.

NEW FEATURES:

- None.

DEFECT REPAIRS:

- Removed the user authentication functionality. A recent update of the authentication platform is causing the authentication to fail preventing our user to complete the installation of the Beeline software successfully.
- Fixed an issue in the Beeline Charts. They now are consistent with the GenomeStudio Controls Dashboard charts.

KNOWN ISSUES:

- None