

Illumina TruSight™ Software Suite Release Notes

v.2.0.0.1

January 2021

Release v.2.0.0.1 Updates

RESOLVED ISSUES

Defect repairs (bug fixes) from v2.0.0 release:

- When a filter that is part of a test definition uses both 'consequence' and 'include genes from phenotype' or saved gene lists, the resulting filter tab will fail for cases using that test definition. Ad hoc filters created after a case is ready for interpretation are not affected.
- While creating a filter during the configuration of a test definition, if the user adds then subsequently removes the gene list sub-filter before saving the filter, the resulting filter tab will not reflect the saved filter for cases using that test definition. There are also certain conditions where the 'include genes from phenotype' filter disappears from the UI after saving the filter. Ad hoc filters created after a case is ready for interpretation are not affected.

KNOWN LIMITATIONS

- Bug fixes above apply to any new filters created in Test Management. Existing filters with these issues will need to be recreated in order to be fixed.

Open known issues from previous releases can be found in TruSight Software Suite Release Notes v.2.0.0 and earlier:

<https://support.illumina.com/downloads/trusight-software-suite-release-notes.html>