

Illumina TruSight™ Software Suite Release Notes

v.2.0.2

August 2021

Release v.2.0.2 Updates

This release includes the fix to resolve intermittent connectivity issues that caused sporadic cases to get stuck in a processing state during ingestion. Users will need to resubmit affected cases to reinitiate ingestion and resume the case processing.

VERSIONS

Versions of component software

- TruSight Software Suite v.2.0.2 runs multiple DRAGEN versions for secondary analysis:
 - Whole genome v.3.5.7b
 - Whole genome v.3.7.5
 - Whole exome v.3.7.5
- TSS v.2.0.2 uses Nirvana v.3.12 for annotation
- TSS v.2.0.2 uses KNS-API v.0.12.0 and KNS-UI v.0.2.0- for curation
- TSS v.2.0.2 uses CLI v.2.0.2

TSS v.2.0.2 uses Emedgene pipeline v.03.55 / 00.05.00.04

RESOLVED ISSUES

Defect repairs (bug fixes) from v2.0.1 release:

- Resolve intermittent connectivity issues that caused sporadic cases to get stuck in a processing state during ingestion.

KNOWN LIMITATIONS

Open known issues from previous releases can be found in TruSight Software Suite Release Notes v.2.0.1 and earlier:

<https://support.illumina.com/downloads/trusight-software-suite-release-notes.html>