

BeadArray Controls Reporter v1.1

Release Notes

Revision History

Revision	Date	Description of Change
v1.1	February 2016	Relaxed the samplesheet constraints & added compatibility with manifest file v1.0+
v1.0	October 2015	Initial Release

Introduction

These Release Notes detail updates following the initial release of the BeadArray Controls Reporter Software. For the software user guide, please see the Documentation & Literature section of the MethylationEPIC array support page. For software download, please see the Downloads section of the MethylationEPIC array support page.

I. BeadArray Controls Reporter v1.1

NEW FEATURES:

- None

IMPROVEMENTS:

- Relaxed the samplesheet constraints. Now Satrix_ID and Satrix_Position are the only mandatory columns.
- Updated BACR formulas for compatibility with Infinium MethylationEPIC manifest file v1.0 B1 and higher.

DEFECT REPAIRS:

- None

KNOWN ISSUES:

- Will not work with Infinium MethylationEPIC manifest file v1.0.